Electric Digital Endodontic Torque Control System

User's Manual

**Manufacturer and Distributor****ATS Dental**
Advanced Technology SolutionsVia Vecchia Provinciale Lucchese, 49
50030 SERRAVALLE PISTOIESE
ItalyIEC 60601-1
EN ISO 13485

This equipment is manufactured according to current regulations and standard (IEC 60601-1) and according to the EN ISO 13485 quality control certification system.

**Service hotline:**

+39 0573 518137

info@atsdental.it

Please indicate the product serial number in all requests. Additional information can be obtained at: www.atsdental.it

This **Operating and Maintenance Instruction Manual** provides instructions on HI TECH ENDO. The text, graphics and icons are used to facilitate understanding of the functions and operations by the final user.

Target Group: Dentists and office personnel.

It is very important to keep this manual with the device, for easy reference at any time in the future. In the event of sale or transfer to another final user, be sure to attach this Operating and Maintenance Instruction Manual. Transferring this manual will ensure, the new final user will have all the necessary information to understand the functionality and proper use of the unit. These instructions are intended to ensure the safe use of this equipment and should be read carefully before installation and use of the device.

Description of safety instruction

	Important information, explanation or integration for users.
	Have a direct impact on the functionality of this product or warns of possible misuse and/or harm to the product or the user.

Three hazard levels are used in this document for avoiding personal and property damage:

 CAUTION	Indicates a hazardous situation that can cause damage to property, or mild or moderate physical harm. (*)
 WARNING	Indicates a hazardous situation that can cause death or serious injury. (*)
 DANGER	Indicates the maximum hazard level. Indicates a directly hazardous situation that can cause death or serious injury. (*)

(*) The optional step covers necessary measures for avoiding hazards.

Symbols

	UE 2017/745		Risk of electric shock
	Follow Instructions for use		Disposal Information
	Type B appliance	V	Supply Voltage
	Manufacturer	REF	Reference
	No access for persons with active implant as pacemaker		Mode: continuous operation with intermittent load
	Safety: Class II		Input Polarity

The symbols printed on the outside are for transportation and storage, and have the following meaning:

	Protect from moisture		Transport upright with the arrows pointing upwards
	Fragile - Protect against impact!	Handle With Care	=

Please keep the packaging in case you need to return the product for servicing or repair.

DFU Rev3 10/2019
HTE800
Pag. 5 / 44

**Do not remove this rating plate from the unit for any reason.
Refer to the SN=Serial Number in any communication with the Manufacturer.**

This symbol placed on the above sticker remind to consult the present
Operating and Maintenance Instruction Manual

I - INTRODUCTION

Congratulations, on purchasing the **HI TECH ENDO** Electric Digital Endodontic System. **HI TECH ENDO** is equipped with the new **Forward - Reverse (Reciprocating) movement** designed for use with files from all major suppliers. Standard **ROTATING** movement is also available.

HI TECH ENDO is a unit with high precision in angle measurement evaluating precisely the angular position.

HI TECH ENDO is very easy to install, operate and maintain.

To receive maximum performance and durability from your **HI TECH ENDO**, please read the entire Operation and Maintenance User Manual before using the device

Separate manuals may be provided to cover the operation and maintenance of handpieces or other accessories for your unit.

Read this instruction booklet carefully before installing and using the machine. In this way, you will obtain the best possible results and maximum operating safety. This booklet must be kept with the apparatus in case of sale or transfer to another user.

II – WARNINGS

To reduce the risks of accidents, the precautions stipulated below must be taken.

Device users:

- Use of the **HI TECH ENDO** is restricted solely to qualified, trained and competent dental health practitioners in the normal context of their work.
- If you have received this device by error, please contact the supplier so that it can be removed.

Interactions:

- The device complies with applicable electromagnetic compatibility standards. The user should nevertheless ensure that any potential electromagnetic interference does not cause an additional risk (presence of radiofrequency emitters, electronic devices, etc.).

Interference may occur when used on patients with cardiac pacemakers. This system emits electromagnetic fields, which means there are some potential risks. The malfunctioning of implantable devices such as cardiac pacemakers and ICDs (implantable cardioverter defibrillator) is possible:

- Ask patients and users if they have an implanted device before using this product.
- Explain the circumstances to them

- Weigh the risks and benefits and contact your patient's cardiologist or appropriate qualified healthcare professional prior to performing the treatment
- Keep this product away from implanted devices
- Make appropriate emergency provisions and take immediate action if patients become ill
- Symptom including a raised heartbeat, irregular pulse, and dizziness may signal problems with a cardiac pacemaker or ICD.

Electrical connection:

- Your device must be connected to the electric power supply by a certified dental installation technician.
- The electric supply to which the device is connected must comply with the standards in force in your country.

Using the device:

- Do not use the device if it appears to be damaged or faulty.
- Turn the device off before unplugging the power cord.
- To unplug the power cord, grip the cord plug and hold the wall socket.
- When the device is not to be used for a long period of time, unplug the device from the electric supply.
- Do not move the device during use.

 CAUTION	<ul style="list-style-type: none">▪ Do not insert anything inside the motor▪ The "Motor" is supplied Non-Sterile!▪ Use only power adaptor provided (chapter 3)
 DANGER	<ul style="list-style-type: none">▪ Use for intended purposes only. Failure to observe the operating instructions may result in the patient or user suffering serious injury. Before using this product, make sure that you have studied and understood the operating instructions.▪ Do not install where there is a risk of an explosion. The systems are not intended for operation in the presence of flammable anesthetics or gases.▪ Do not disassemble or alter the system motor, console, or footswitch.
 WARNING	<ul style="list-style-type: none">▪ Connect mains cord to a properly outlet only.▪ Never touch files when they are still rotating.▪ Contra-angle should only be attached when the motor has stopped running.

Environment:

- Do not cover the device.
- Do not immerse the device in liquid, and do not use it outdoors.
- Do not tilt the device at an angle greater than 5°.
- Do not place the device near a heat source.
- Do not change leads with leads different than the supplied ones.
- Make sure that the cords are not in a traffic path.

- The device should be stored in its original packaging in a safe place.
- The device is not designed for use in the presence of anesthetic gases or any other flammable gas.
- Do not expose the device to water vapor, or splashes.
- The device is not designed to work near ionizing radiation.
- Do not insert metal objects into the device (risk of electric shock, short-circuit or emission of hazardous substances).

Maintenance:

- Before and after each use, your device must be disinfected with products recommended by the manufacturer.
- Before each procedure, it is essential to make sure that the accessories to be used have been cleaned and disinfected.

Accessories:

- The device can only accept accessories designed by the manufacturer for the particular use for which they are intended.
- Use of accessories from other manufacturers is a potential hazard for you or your patients.

Repair:

- Do not repair or modify the device without prior authorization from the manufacturer.
- In the case of a fault, contact the supplier of your device. Do not use unauthorized repairers, who might make your device dangerous for you and your patients.


If you have any doubt, contact an approved dealer or our customer support department:

- E-mail: info@atsdental.it


III – DESCRIPTION

PHYSICAL DESCRIPTION

1. HT801 - Electronic Control Console
2. HT802- Micro Motor
3. HT803 - Foot Control
4. HT804 - Power Supply9V 2A

TECHNICAL DESCRIPTION

**TFT screen/control Keypad**

The HI TECH ENDO is adjusted by applying reasonable pressure to the keypad buttons.

1. "MODE" The operator can easily switch from "F+R" and "ROTATING".
2. "PROGRAM" The Display change accordingly. Scrolls Programs (1 -10).
3. "RATIO" Set the handpiece ratio. (1:1 - 6:1 - 16:1 - 20: 1- 32:1.)
4. "CALIBR." Performs a self-calibration of handpiece and micromotor.
5. "MOTOR" Switches ON and OFF motor (if enabled).
6. "SPEED+/-" Set motor rotation speed (value in rpm).
7. "TORQUE+/-" Set motor torque (value in Ncm).
8. "FWD/REV" Clockwise and anti-clockwise.
9. "REV type" Stop, Auto reverse, Special Auto Reverse.
10. "FWD°" Scroll up and down from 10° to 360 ° in forward.
11. "REV°" Scroll up and down from 10° to 360 ° in Reverse.

Information required by the user is displayed on the TFT screen.

b. Control unit back panel

1. Main Power On/Off Switch.
2. Power

c. Control footswitch

Activates the pre-selected maximum speed or start Forward/Reverse Movement

d. Technical characteristics

Console / micro motor / footswitch			
Manufacturer	ATS DENTALsrl	Fuses	NO
Dimensions (W x H x D)	170 x 110 x h90 mm	Protection	Class II
Materials	ABS	Type	B
Micro motor	10000 rpm / Torque: 6.8 mNm	Class MD	Ila rule 9 in Annex VIII UE2017/745
Noise	< 65 dBA	Protective index console	IP20
Input Voltage	9 VDC	Protective index micro motor	IP65
Input Current	2.0 A	Protective index footswitch	IP20
Not suitable for use in presence of flammable anesthetics or oxygen.			
Conditions of use: Temp +18°C/+40°C (+64°F/+140°F) RH < 80%. Altitude restriction:3.000m (10.000 ft)			
Shipping and Storage conditions: Temp +5°C/+65°C (+41°F/+149°F) RH < 20-95% non condensing.			
Life of the device: When used as prescribed = 5 years. No warranty claim may be made in the event that the failure happen sooner or later than indicated, because of the frequency of use, sterilization and maintenance.			

IV– INSTALLATION / FIRST USE**UNPACKAGING THE DEVICE**

Upon reception of the device, check for any damage caused in transit.

Contact your supplier if necessary.

RECOMMENDATIONS

Have your device connected to the power supply by a certified dental installation technician. The electrical connection of the HI TECH ENDO must comply with the standards applicable in your country.

INSTALLATION**Important:**

Do not place the HI TECH ENDO close to or on top of another device.

Do not place the power cord and the footswitch cord in a wire cover or in a cable gland.

1. Remove the unit and its accessories from their original packaging and place on a flat surface. Do not install near heat sources, direct or indirect.
2. Verify the contents of the box according to the section “**PHYSICAL DESCRIPTION**”.

3. Do not kink the motor hose since it may damage it.
4. All accessories are supplied non-sterile.
5. The serial number on the unit's rear panel must be the same as the one on transport documents.
6. Allow at least 150 mm free space around the unit for cooling ventilation.
 - Place the control unit in position on a fixed surface, horizontal or with a slope of no more than 5°.
 - Check that the power switch is in position O (off).
 - Connect the power cord to the receptacle of the device.
 - Connect the power cord to a socket.
 - Connect the control footswitch to the footswitch connector.
 - Position the footswitch so that it is easily accessible for your feet.
 - Connect the micro-motor and cord to the connector.
 - Place the micro-motor on its rest.

V–ADJUSTMENTS / SETTINGS / INTERFACE

START UP

Before turning on the central unit, check to see that you have completed the installation of ALL COMPONENTS.

When first switched on the unit will show “**Rotating**”, and values set by manufacturer on the display.

Use “MODE” button to toggle from “Rotating” Mode to “F+R” Mode

SETTINGS

HI TECH ENDO will store the operative parameters for each “program” as soon as they will be set by the user and they will be maintained in the memory even if the central unit will be switched off. For each “program” it is possible to set: speed, torque, ratio, type of auto-reverse, values of the Angle in FWD and REV. **To reset the parameters to the manufacturer ones** (Software unit reset), **press the “Program” key while the unit is switched on.**

TORQUE CALIBRATION

Perform torque calibration, press “**Calibration**” key on front panel while holding the motor. This function check the handpiece mechanical inertia, in order to off-set this value. This operation will take a few seconds, the motor will run at different speed (Fig. 7)

Time of Calibration: about 8 sec. Do not perform any operation during calibration. Calibration should be done before each endo session or anytime the contra-angle is changed or lubricated.

 CAUTION	<ul style="list-style-type: none">The file must not be engaged in any operation when this function is performed.Calibration should be done before each session or anytime the handpiece is changed or lubricated.
 i	<ul style="list-style-type: none"><i>Do not perform calibration if the micro motor is not inserted in the receptacle on the front panel of the unit. In case of error, perform a "System Reset" as described.</i>

PROGRAM SETTINGS ADJUSTMENT

Rotating mode

In every program has been set the same information by the manufacturer:

The listed parameters are only examples. To prevent unnecessary risk, observe the guideline speeds given by the manufacturer of the rotating instruments and contra-angle system in use.

Programming the unit:

- RATIO:** Sets handpiece reduction ratio. See the available reduction ratio in the following table:

1:1	6:1	16:1	20:1	32:1
-----	-----	------	------	------

2. SPEED +/-: Sets motor rotation speed. See the following table:

	Reduction Rate				
	1:1	6:1	16:1	20:1	32:1
RPM	10000	1660	600	500	300
	9000	1500	550	450	280
	8000	1330	500	400	250
	7000	1160	450	350	200
	6000	1000	400	300	180
	5000	800	300	250	150
	4000	650	250	200	100
	3000	500	200	150	
	2000	300	100	100	
	1000	150			

3. TORQUE +/-: increases or decreases the torque values (24 steps). Torque values are in Ncm and depend on contra-angle reduction rate. See the following table:

Ncm with Ratio 1:1 or 6:1						
0.5	1.0	1.5	2.0	2.5	3.0	3.5

Ncm with Ratio 16:1 20:1 32:1						
0.5	1.0	1.5	2.0	2.5	3.0	7.5

4. REVERSE: Switch between forward and reverse motor rotation (when in reverse, console will emit a beeping tone).

5. REV TYPE :Set the type of Auto-reverse:

A. STP – Auto-reverse “OFF”:

The motor handpiece stops without reverse rotation

B. AR – Autoreverse “Standard”:

When file reaches the set torque value, the motor runs reverse until the motor has been stopped. Next start will be in Forward. *NB: In use only if Torque is not set on Max*

C. SAR – Auto-reverse “Special”:

When file reaches the set torque value, the motor runs reverse until it has freed the file.
Then the motor runs again in forward. *NB: In use only if Torque is not set on Max*

NB:

Beeping sound notifies user when 75% of the preset torque has been reached and when preset torque limit has been exceeded.

CAUTION

Always check the values before use. Danger of injury for incorrect use.

F+R mode

The listed parameters are only examples. To prevent unnecessary risk, observe the guideline speeds given by the manufacturer of the rotating instruments and the guideline given by the manufacturer of the contra-angle in use.

Programming the unit:

1. **RATIO:** Sets handpiece reduction ratio. See the available reduction ratio in the following table:

1:1	6:1	16:1	20:1	32:1
-----	-----	------	------	------

2. **RPM +/- :**Sets motor rotation speed. See the following table:

	Reduction Rate				
	1:1	6:1	16:1	20:1	32:1
10000	1660	600	500	300	
9000	1500	550	450	280	
8000	1330	500	400	250	

RPM	7000	1160	450	350	200
	6000	1000	400	300	180
	5000	800	300	250	150
	4000	650	250	200	100
	3000	500	200	150	
	2000	300	100	100	
	1000	150			

3. **FWD° +/-:** increases or decreases the angle value in **FWD** from 10° to 360°, by step of 5°.
4. **REV° +/-:** increases or decreases the angle value in **REV** from 10° to 360°, by step of 5°.

Always check the values before use. Danger of injury for incorrect use.

VI – USING THE DEVICE

Before turning on the central unit, check to see that you have completed the installation of ALL COMPONENTS.

Important:

- Do not disconnect the micro-motor cord when the device is switched on and the footswitch is pressed.
- Do not engage the instrument to or disengage it from the contra-angle when the micromotor is running.
- All accessories must be cleaned and disinfected before use.
- For your safety and that of your patient, the HI TECH ENDO must not be used with accessories other than those supplied by the manufacturer.
- Check the integrity of the device and its accessories before and after each use in order to detect any problem.
- If necessary, do not use the device and replace any faulty component.

The HI TECH ENDO is prepared for use as follows:

- Make sure that the device is connected correctly and properly insulated.

- Micromotor connection:

- Connect the micro motor into the receptacle on the front of the central unit
- Connect the contra-angle to the micro motor and fit an instrument (file) to the contra-angle.

Attach the appropriate "E-Type" handpiece to the motor. It can be used with most endodontic contra-angles without internal light or spray on the market.

- Use only contra-angle that conform to ISO 3964.
- Change the contra-angle only when the motor not running.

- Place the micro motor on the motor holder.

- Footswitch connection:

- Connect the footswitch into the receptacle on the front of the central unit
Place the footswitch on the ground.

While inserting the plug, align the guide on the connector in an upright position to fit into the corresponding female guide on the HI TECH ENDO (Fig. 13). Connectors have a spring retention system to help avoid accidentally disconnecting the unit's cords and cables.

- Switch on the device:

Electrical connections

- Check the voltage present on your power line (100VAC to 240VAC).
- Connect the power cord to HI TECH ENDO by inserting the female connector of the power supply into the male receptacle on the back of the unit.
- Connect the power supply into the electrical outlet.

VII – SHUTTING DOWN THE DEVICE

Upon completion of the dental procedure:

- Remove the instrument fastened to the contra-angle.
- Disconnect the contra-angle from the HI TECH ENDO micro-motor.
- Disconnect the micro-motor cord from the HI TECH ENDO unit.
- Switch the device off (O).

VIII – ROUTINE MAINTENANCE / STERILIZATION

Important:

The HI TECH ENDO micro-motor and cord is not autoclavable.

CONTRA-ANGLE MAINTENANCE

Refer to the manufacturer's instructions for your contra-angle.

ROTARY INSTRUMENT MAINTENANCE

Please refer to the manufacturer's instructions for your instruments (file).

DEVICE MAINTENANCE

Important:

Do not use an abrasive product to clean the device.

Do not allow cleaning/disinfecting spray to penetrate inside the device.

The device, motor and footswitch are not sterilizable.

Unplug the main unit from power source before cleaning or disinfecting. Do not use directly on liquids and sprays. Never immerse the unit or its parts in liquid.

External surface of the console may be cleaned with a damp cloth or disinfectant wipe. The housing is not waterproof. Do not use ultrasonic or steam sterilizers.

The HI TECH ENDO control unit and the footswitch must always be cleaned and disinfected after each procedure using alcohol, disinfectant agents or disinfectant wipes for dental practices.

The exterior of the footswitch may be cleaned by wiping with a soft cloth moistened with mild detergent or disinfecting solution.

Regular monitoring of the HI TECH ENDO unit is necessary in order to detect any problem. It is important to keep the control unit ventilation vents clean in order to avoid abnormal heating.

MICRO-MOTOR AND CORD MAINTENANCE

Do not sterilize or lubricate the micro motor.

Do not attempt to disassemble the motor or motor connector. Unplug the main unit from power source before cleaning or disinfecting. Do not use directly on liquids and sprays. Never immerse the motor in liquid. External surface of the micro motor may be cleaned with a damp cloth or disinfectant wipe. Do not place micro motor in ultrasonic sterilizers. Motor holder may be cleaned with a damp cloth or disinfectant wipe.

 CAUTION	Always check the values before use. Danger of injury for incorrect use.
	<i>When treating patients who may have an acute, critical infectious disease, be sure to observe the hygienic measures cited in applicable publications and reports. If possible, use suitable disposable products to avoid the transmission of critical pathogens. These protect the user, the patient and all participants in the office.</i>

IX – MONITORING / PREVENTIVE AND CORRECTIVE MAINTENANCE

Important:

In the case of a fault, you are advised to contact the supplier of your device. Do not use unauthorized repairers, who might make your device dangerous for you and your patients.

MONITORING

Regular monitoring of the device and its accessories is necessary in order to detect any insulation fault or damage. Replace, if necessary.

It is important to keep the control unit ventilation vents clean in order to avoid abnormal heating.

PREVENTIVE AND CORRECTIVE MAINTENANCE

HI TECH ENDO micro-motor inspection/overhaul: The manufacturer recommends inspection or overhaul of the I-Endo micro-motor at least once a year.

The micro motor cable has limited durability and is considered a consumable product. The durability of the cable is largely dependent upon use and treatment, i.e., handling, reprocessing and frequency of use.

Caution must be exercised when connecting or disconnecting the micro motor and foot control cable from the device to avoid damage to the cable (Fig. 13).

Periodically inspect the micro motor cord for damage prior to use. If there is damage noted, contact ATS DENTAL.

1. Do not attempt to disassemble the motor or motor connector.
2. Do not oil or lubricate the motor.
3. Do not attach a handpiece to the motor while the motor is running.
4. Do not bend motor cord sharply.
 - The micro motor is sensitive to shock. Do not drop or impact micro motor against hard surface

Failure to comply with any of the above instructions may void your warranty

The average life of this medical device with normal use and in accord with the present user manual, has been calculated to be 5 years. According with EN62353 we recommend users verify every 2 years from the date of the first use: dispersion current versus ground, in the box and in the patient.

MOTOR RUBBER O-RING:

- If necessary, substitute the "Rubber O-Ring"

LUBRICATING:

- **Do not oil or lubricate the micro motor.**

	<i>Micro motor and console do not require modification.</i>
---	---

OPERATION FAULTS**Fault detected: The central unit does not switch on.**

Possible cause	Solution
• Power supply not connected properly	Make sure the power supply is properly plugged into the electrical outlet and the connector on the back of the unit.
• Unit has not been switched on.	Check the back of the unit to see that the On button has been pressed.
• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another working electrical outlet. If the unit does not turn on, verify the integrity of the fuse in
• Broken connections inside. • PCB faulty. • Power interrupted.	Immediately discontinue use of the device and contact ATS DENTAL after-sales-service.

Fault detected: The central unit turns off and stops working suddenly

Possible cause	Solution
• Power supply cable damaged. • PCB damaged. • Connections inside interrupted	After checking the correct connection to the electrical outlet and the absence of liquid in the unit, try to power the central unit. If the device does not turn on, immediately discontinue use of the device and contact ATS DENTAL

Fault detected: The foot pedal doesn't start the motor

Possible cause	Solution
• Micro Motor cable not connected properly	Check to see that the micro motor cable is securely plugged into the connector on the front right of the central unit, disconnect and reconnect.
• Footswitch cable is not connected properly	Check that the footswitch cord is plugged into the connector on the front left of the central unit, disconnect and reconnect
• Male connector of the micro motor damaged.	Unplug the micro motor and make sure that the male connector is not damaged, wet or damp.

• Pins of the connector damaged	Unplug the micro motor and the footswitch cable, make sure the metal pins on the male connector are intact and not damaged and have not broken off and/or remained stuck in the female connector mounted on the unit.
---------------------------------	---

• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
------------------	---

• Footswitch damaged.	Check the ability to start the micro motor from the central unit keyboard, pressing the [Motor ON/OFF switch]. If so, contact ATS DENTAL after-sales-service.
-----------------------	---

Fault detected: The display is not working properly

Possible cause	Solution
• Contra angle damaged.	Check the integrity of the contra angle.
• Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
• Possible interference.	Check to see if the central unit has been positioned too close to another electrical device which can cause interference. In this case, move the interfering electrical device further away from the central unit.
• Broken connections inside. • Display damaged. • PCB damaged.	If you cannot find the cause of the malfunction, immediately discontinue use of the device and contact ATS DENTAL after-sales-service.

Fault detected: The keyboard is not working properly

Possible cause	Solution
• Keyboard damaged.	Contact ATS DENTAL after-sales-service.
• Broken connections inside.	

Fault detected: The sound signals are not working properly

Possible cause	Solution
• Broken connections inside	Immediately discontinue use of the device and contact ATS DENTAL after-sales-service.

Fault detected: The motor doesn't work

Possible cause	Solution
• Micro motor cable not connected properly	Check to see that the micro motor cable is plugged into the connector on the front right of the central unit, disconnect and reconnect.
• Contra angle defective or not properly calibrated.	Check the integrity of the contra angle. Remove the contra angle and perform calibration.

•Footswitch damaged.	Attempt to start the micro motor by using the [Motor ON/OFF button] on the central unit keyboard. If unit does not turn on, contact ATS DENTAL after-sales-service.
•Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
•PCB damaged. •Motor damaged. •Connector pins and/or cable damaged.	Immediately discontinue use of the device and contact ATS DENTAL after-sales-service.

Fault detected: The motor stops suddenly

Possible cause	Solution
•Micro motor/Cable damaged	Check to see that the micro motor cable is plugged into the connector on the front right of the central unit, disconnect and reconnect. Start the motor by pressing the [Motor ON/OFF button on the keyboard] . If the motor turn on, immediately discontinue use of the device and contact ATS DENTAL after-sales-service.
•Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.
•Footswitch damaged.	Attempt to start the micro motor by using the [Motor ON/OFF button] on the central unit keyboard. If unit does not turn on, contact ATS DENTAL after-sales-service.
•PCB damaged.	If you cannot find the cause of the malfunction, immediately discontinue use of the device and contact ATS DENTAL after-sales-service.

Fault detected: The Motor will not turn off

Possible cause	Solution
•Keyboard and/or footswitch damaged	Try to turn off the motor by pressing the [Motor ON/OFF button] on the keyboard or by using the footswitch. If the motor will not turn off, immediately discontinue use of the device and contact ATS DENTAL after-sales-service.
•Broken connections inside. •PCB damaged	Immediately discontinue use of the device and contact ATS DENTAL after-sales-service.
•Wrong voltage.	Verify that the voltage matches the power supply. If it does, try to switch the central unit off and then on again, or connect the device to another electrical outlet.

X – ELECTROMAGNETIC COMPATIBILITY

Important:

The power cord, micro-motor cord, and footswitch cord must be kept apart.

The HI TECH ENDO requires special precautions to be taken with regard to electromagnetic compatibility. It must be installed and prepared for use as described in chapter 4.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the HI TECH ENDO. The recommended separation distances in this paragraph must therefore be complied with.

The HI TECH ENDO must not be used near to or on top of another device. If this cannot be avoided, its operation under the conditions of use must be checked beforehand.

The use of accessories other than those specified or sold by the manufacturer as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the HI TECH ENDO.

 CAUTION	To avoid possible risks of electromagnetic interferences do not use other electro-medical instruments of similar nature near HI TECH ENDO. The unit complies with the current normative electromagnetic radiation law.
	This device has been tested and found to comply with the emissions requirements of IEC 60601-1-2. These requirements provide reasonable protection against harmful electromagnetic interference in a typical medical installation. However, high levels of radio-frequency (RF) emissions from electrical devices, such as cellular phones, may disrupt the performance of this device. To mitigate disruptive electromagnetic interference, position this device away from RF transmitters and other sources of electromagnetic energy.
 DANGER	Risks from electromagnetic fields: The functions of implanted systems (such as pacemakers) can be influenced by electromagnetic fields. Patients and users must inquire about implanted systems before using themotors and check the use. Do not bring the motors close to the systems. Take suitable emergency precautions and immediately react to any changes in health.

ELECTROMAGNETIC COMPATIBILITY

HI TECH ENDO is intended for use in the electromagnetic environment specified below. The HI TECH ENDO should only be used in such an environment.

RF emissions CISPR 11	Group 1	HI TECH ENDO uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	HI TECH ENDO is suitable for use in all establishment including domestic establishments and those directly connected to the public low-voltage power supply.
Harmonic emissions IEC 61000-3-2	Class A	
Immunity ESD	15kV in air 8kV contact	Floor should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%
burst	2kV/100kHz	
Magnetic field:	30A/m	
immunity to RF currents	range 150kHz- 80MHz	3V modulation 80% 1kHz
surge immunity		1kV differential mode 2kV common mode

XI – DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 23/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the nearest ACTEON GROUP office, the list of which is given in chapter 15), for information on how to proceed.

XII – MANUFACTURER'S LIABILITY

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.),

- repairs or other work have been done on the device by persons not authorized by the manufacturer,
- the device has been used in an electrical installation which does not comply with current standards,
- the device has been used in a way which is not stipulated in this manual,
- accessories (motor and cord, footswitch, etc.) other than those supplied by ATS DENTAL have been used,
- the instructions in this document have not been followed.

Note:

The manufacturer reserves the right to modify the device and/or the user manual without notice.

XIII – REGULATIONS

This medical device is classified as class IIa according to European directive 93/42/EEC. This equipment is manufactured in compliance with the current IEC 60601-1 standard. This equipment has been designed and manufactured according to an ISO 13485-certified quality assurance system.

ATS riserve the right to change this manual without advice

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