


THE NEW GENERATION....OF POWER AND CONTROL

I-TECH Lite

Code «ITL500»

Surgical Micromotor System Touch Screen - Safe - Powerful



Autoclavable Brushless Micromotor
Guaranteed up to 1000 cycles




40000 RPM



Torque: 80 Ncm



135°C
Light without shadows
Power: 31000Lux



Electronic High Precision Calibration
Ensures safe treatment
trought accurate torque correction




**LED
Lite**

Wide light area
Light without shadows
Power: 31000Lux

«DeLuxe» Multifunction Foot Control

Covered entirely with syntetic rubber: resistant to any leaking fluid. The reostatic effect is ensured by means of a specialized «pressure sensor». There are NO mechanical or moving parts.

4.3" TOUCH SCREEN



Advanced Irrigation Pump

Easy set up of irragation tubes - Very ergonomic opening and closing system - Better control of irrigation flow with very constant flow rate - Stable and quiet operations - 3 Flow Rates available: MIN - MED - MAX



INCLUDED

