

«Omega» Implant



Brushless Micromotor



Electronic Calibration of ContraAngle



Available with 3 different Pedal



Built to Last



User Friendly




30000 RPM



Torque: 51 Ncm

**CA
20:1**




**The new generation...
of Power and Control**

**Available with
Constrangle 20:1**



«Omega» Implant

Implant Motor with integrated
coolant pump easy and uncomplicated to use
Brushless motor



Our **OMEGA IMPLANT** is a compact tabletop unit. All necessary treatment stages of implantology can be performed safely and precisely using just one 20:1 contra-angle.

OMEGA IMPLANT is equipped with **Contra-Angle calibration**: the precision of the motor is adapted to the contra-angle connected. **OMEGA IMPLANT** also delivers the highest standards of ergonomics and flexibility.

The motor torque is **72 Ncm** and the motor speed range is **400 – 40,000 rpm (20 - 2000 with 20:1 contra-angle)**.




«DeLuxe» Multifunction Foot Control

Programs switching - Activation and Deactivation of the coolant pump - Clockwise/ Anticlockwise rotation - ON/OFF function. The reostatic effect is ensured by means of a **specialized «pressure sensor»**. There are NO mechanical or moving parts.



«Standard» Multifunction Foot Control

Programs switching - Activation and Deactivation of the coolant pump - Clockwise/ Anticlockwise rotation - ON/OFF function. with **«ERC» (Electronic Reostatic Effect)**



«Easy» Foot Control

ON/OFF button with **«ERC» (Electronic Reostatic Effect)**

Advanced Irrigation Pump



- Easy set up of irrigation tubes - Very ergonomic opening and closing system - Better control of irrigation flow with very constant flow rate - Stable and quiet operations - 3 Flow Rates available: MIN - MED - MAX

ATS Dental


Via Vecchia Provinciale Lucchese, 49
51034 SERRAVALLE PISTOIESE (PT)-ITALY
info@atsdental.it - www.atsdental.it

Model/Code	Description
AOMG300DL	Implantology Unit:
AOMG300S	- Electronic Control Console
AOMG300E	- Brushless Micromotor
	- Foot Control
	- Solution Holder
OMEGA IMPLANT	- Motor Holder - 2 Tube