

«Implant» Accessories

● Without Light



Autoclavable, Built to Last, User Friendly

● With Light




Autoclavable, Built to Last, User Friendly 20:1 Contra-Angle - Push Button - Fiber Optic. Excellent light condition facilitate accurate diagnoses and perfect treatment results



	CODE: 201-ATS	CODE: 201-ATSL
<i>Indications</i>	Implantation Osteotomy	Implantation Osteotomy
<i>Transmission Ratio</i>	20:1	20:1
<i>Coupling System</i>	ISO 3954/DIN 13940	ISO 3954/DIN 13940
<i>External Irrigation</i>	✓	✓
<i>Internal Irrigation</i>	✓	✓
<i>Chucking System</i>	push-button	push-button
<i>Lighting LED System</i>	no	✓
<i>Rotary Instruments</i>	for surgical burs shank Ø 2.35 mm	for surgical burs shank Ø 2.35 mm
<i>Max Drive Speed</i>	40,000 min ⁻¹	40,000 min ⁻¹

ATS Dental



Via Vecchia Provinciale Lucchese, 49
51034SERRAVALLE PISTOIESE (PT) - ITALY
info@atsdental.it - www.atsdental.it


Code	Description
201-ATS	20:1 Contra Angle without Light
201-ATSL	20:1 Contra Angle with Light



«Implant» Accessories

IRRIGATION LINE ...Surgical irrigation tubing set, single-use, sterile.

Sterile, flexible, external tubing **with integrated peristaltic pump**, intended to provide a conduit to deliver a sterile irrigation solution from its source to the surgical intervention site.



Code: ATSDEF1

Internal or external mechanical irrigation Kit for surgical physiodispenser

CONF: 10 pz

Code	Description
ATSDEF1	Irrigation Line 1 Via
ATSDEF2	Irrigation Line 2 Vie

CARRYING CASE



Code: ATSCC

The ATS Carrying Case can accommodate all «Implant Units» components

Code	Description
ATSCC	Carrying Case



Via Vecchia Provinciale Lucchese, 49
51034SERRAVALLE PISTOIESE (PT) - ITALY
info@atsdental.it - www.atsdental.it

