

# LAB X20

## HEAVY DUTY LABORATORY MOTOR

- For heavy acrylic use
- Powerful one piece handpiece
- 35k rpm micromotor
- Console with adjustable speed and direction

**LAB X20** Electric Laboratory Handpiece is a compact electric lab handpiece system that provides excellent cutting power and fully variable speed from 0 to 35,000 rpm.

Designed for all dental lab applications such as cutting, grinding, sanding, trimming, etching, polishing and more, the **LAB X20** features a high torque and a lightweight motor. The motor allows reverse direction and a cruise function provides consistent rotation speed.

Model	Code	Description
LAB-X20	LBX20	Heavy Duty Lab Motor: - Electronic Control Console - Handpiece - Foot Control



**ATS** Dental



Via Vecchia Provinciale Lucchese, 49  
51034SERRAVALLE PISTOIESE (PT) - ITALY  
info@atsdental.it - www.atsdental.it

